

CS-200M

Chair Scale



- Capacity– 200 kg / 440 lb
- Dual Marked in Kilograms / Pounds
- Locking Rear Wheels for Safety
- Durable Construction
- Swivel Front Wheels for Easy Transport
- Accurate Full Range Weights
- Easy to Read Dual Beam
- Lift Away Armrests for Easy Access

VALUE IN
WEIGHING

CS-200M Chair Scale

Features

This mechanical chair scale was designed for weighing those who have difficulty standing on their own. The die-cast easy-to-read beam runs along the back of the chair. Powder coated steel construction.

The chair is oversized for comfort with armrests that lift away for easy access. The rear wheels lock for safety and the front wheels swivel for easy portability along with easy-to-clean plastic seat.

Specifications

Capacity: 200 kg x 100 g / 440 lb x 4 oz

Seat Dimension:

405 mm (W) x 350 mm (D) x 340 mm (H)

16" W X 14" D x 13" H

Product code: 816965004843

Product Weight: 32 kg / 70 lb

Shipping Weight: 36 kg / 80 lb

Package Dimension:

762 mm (W) x 609 mm (L) x 635 mm (H)

30" W x 24" L x 25" H

Packaging Contents: Manual, footrest, 2 armrests, scale with wheel assembly, seat, back handle, weighing tower and miscellaneous hardware



a division of Avery Weigh-Tronix, LLC

Brecknell

1000 Armstrong Dr., Fairmont, MN 56031
Toll free: 800-637-0529, Tel: 507-238-8702
Fax: 507-238-8271
e-mail: sales@brecknell scales.com
www.brecknell scales.com



Brecknell

Foundry Lane, Smethwick
West Midlands B66 2LP UK
Tel: +44 (0) 845 246 6717
Fax: +44 (0) 845 246 6718
Email: sales@brecknell scales.co.uk
www.brecknell scales.co.uk

