

HS-200M

Physician Scale



- 200 kg / 440 lb Capacity
- Height Gauge to 210cm / 7'
- Rollers Included for Easy Transport
- Easy to Read Dual Beam
- Dual Marked in Kilograms / Pounds
- Durable Construction
- Accurate Full Range Weights
- Large 375 mm x 275 mm / 14.75" x 10.8" Weighing Platform

VALUE IN
WEIGHING

HS-200M Physician Scale

Features

The HS-200M is a mechanical device for use where either weight or height measurements are necessary. The scale provides accurate weight or height measurements on easy to read linear beams in bold marking. It has a large weighing range and can measure height from 60 cm to 210 cm x 0.1 cm / 2' to 7' x 0.125".

Constructed of steel with a durable white finish, this scale will provide years of service in normal applications.

Specifications

- Capacity:** 200 kg x 100 g / 440 lb x 4 oz
- Height Rod Range:**
60 cm to 210 cm x 0.1 cm / 2' to 7' x 0.125"
- Weighing Platform Size:** 375 mm (L) x 275 mm (W)
14.75" L x 10.80" W
- Overall Dimension:**
- | Length | Width | Height |
|--------|--------|---------|
| 559 mm | 275 mm | 1485 mm |
| 22" | 10.8" | 58.5" |
- Product Code:** 816965004850
- Product Weight:** 17 kg (38 lb)
- Shipping Weight:** 18 kg (40 lb)
- Package Dimension:**
1499 mm (L) x 508 mm (W) x 152 mm (H)
59" L x 20" W x 6" H
- Packaging Contents:** Manual, platform, height scale, weighing tower, wheels and misc hardware



a division of Avery Weigh-Tronix, LLC

Brecknell
1000 Armstrong Dr., Fairmont, MN 56031
Toll free: 800-637-0529, Tel: 507-238-8702
Fax: 507-238-8271
e-mail: sales@brecknellscales.com
www.brecknellscales.com



Brecknell
Foundry Lane, Smethwick
West Midlands B66 2LP UK
Tel: +44 (0) 845 246 6717
Fax: +44 (0) 845 246 6718
Email: sales@brecknellscales.co.uk
www.brecknellscales.co.uk

